

THE BIPOLAR 5 mm OFFICE RESECTOSCOPE



Click or scan here
to watch the teaser video

Versatile use

- Minimal outer diameter
 - Optimized intervention time*
 - Increased safety, accessibility and patient satisfaction*
 - Possible to perform resection in the doctor's office or clinical outpatient setting, even without anesthesia
 - Little or even no dilation required (depending on the patient)
- Optional use as „see & treat“ hysteroscope
- Real bipolar system with current not returned via the sheath

* <http://www.fvvo.be/assets/725/FVWinObGyn-10-115.pdf>.









www.karlstorz.com




STORZ
KARL STORZ – ENDOSKOPE

THE BIPOLAR 5 mm OFFICE RESECTOSCOPE



Due to its small diameter, this resectoscope offers the possibility to perform resections outside the OR. Infertility patients or those suffering from abnormal bleeding can thus be treated in the doctor's office or clinical outpatient setting, even without anesthesia.

Optic and Working Element Set	
26120 AA	HOPKINS® Straight Forward Telescope 0° , diameter 2.9 mm, length 30 cm
26053 SCK	Resectoscope Sheath , 15 Fr., oblique beak, rotating inner sheath, color code: green
26053 OC	Standard Obturator
26053 EBH	Working Element Set , bipolar
26053 CD	Telescope Bridge , with channel for semi rigid 5 Fr. operating instruments, for use with Resectoscope Sheath 26053SC

HF Electrodes for various indications (sterile, for single use, package of 10)			
	011050-10*	Cutting Loop , bipolar, 15 Fr.,	STERILE 
	011051-10*	Cutting Electrode , bipolar, pointed, 15 Fr.,	STERILE 
	011052-10*	Coagulation Electrode , bipolar, ball end, 15 Fr.,	STERILE 
	011053-10*	Cutting Loop , bipolar, longitudinal, 15 Fr.,	STERILE 

Cold Loops for blunt dissection/enucleation or retrieving resected tissue		
	26053R	MAZZON Cold Loop , straight, 15 Fr.
	26053T	MAZZON Cold Loop , rake-shaped, 15 Fr.
	26053U	MAZZON Cold Loop , knife-shaped, 15 Fr.

HF Device	
UH400	AUTOCON® III 400 , high-end HF generator, without CE certification, power supply 220--240 VAC, 50/60 Hz
UH400E	Same , with CE certification

Pump Systems	
26340001-1	HYSTEROMAT E.A.S.I.® SCB
031717-10*	Tubing set , irrigation STERILE 
031217-10*	Tubing set , suction STERILE 



For direct contact email us at gyn@karlstorz.com

It is recommended to check the suitability of the product for the intended procedure prior to use. Please note that the described products in this medium may not be available yet in all countries due to different regulatory requirements.



KARL STORZ SE & Co. KG
Dr.-Karl-Storz-Straße 34, 78532 Tuttlingen/Germany
Postbox 230, 78503 Tuttlingen/Germany
Phone: +49 7461 708-0
Fax: +49 7461 708-105
E-Mail: info@karlstorz.com

www.karlstorz.com

